

**With LED light**



# SCALEX® BRITE


## PIEZO ULTRASONIC SCALER WITH LED

- Autoclavable handpiece with LED light technology.
- Built-in LED allows for direct illumination of the oral cavity.
- Futuristic touch sensitive panel technology.
- State-of-the-art piezoelectric technology is self-tuning and load-adjusting for smooth and effortless operation.
- Advanced water-cooled circuitry designed for greater stability during operation.
- Extremely gentle on the patient.
- Wide selection of high quality tips to correspond with specific clinical applications.
- Available for 100-240V 50/60Hz operation.

**IT'S ALL IN THE HANDS.**

**DENTAMERICA®**

## VALUE BY DESIGN



### PIEZO SCALER TIPS



88002

Reorder No.155-03  
Endo Chuck



88003

Reorder No.155-01  
Universal A



88004

Reorder No.155-51  
Universal B



88005

Reorder No.155-02  
Chisel



88006

Reorder No.155-48  
Perio

### SPECIFICATIONS

INPUT: AC 100V~240V 50/60Hz

OUTPUT: DC 24V 1.5A

FREQUENCY: 27,000 to 32,000 Hz

DIMENSIONS:  
W221mm x D142mm x H43mm

WEIGHT: 620g

DIMENSIONS(HP): Ø18.5mmxW105.5mm

WEIGHT(HP): 49g

ISO 13485 **GMP**

#### DENTAMERICA INC.

18688 E. San Jose Ave. City of Industry, CA 91748 U.S.A.  
TEL: (626) 912-1388 FAX: (626) 913-0510  
www.dentamerica.com  
www.dentamerica.asia

DISTRIBUTOR

**DENTAMERICA®**