



# SCALEX® 980

## PIEZO ULTRASONIC SCALER

- Futuristic touch sensitive panel technology, easy cleaning and maintenance.
- Autoclavable handpiece also applies to endodontic tip.
- State-of-the-art piezoelectric technology is self-tuning and load-adjusting for smooth and effortless operation.
- Advanced water-cooled circuitry designed for greater stability during operation.
- Special TURBO function helps remove especially hard calculi by increasing power relative to the hardness.
- Extremely gentle on the patient.
- Wide selection of high quality tips to correspond with specific clinical applications.
- EMS type tips compatible.

**BEST SCALING WITH A LIGHT TOUCH.**

**DENTAMERICA®**

# VALUE BY DESIGN



Order No. 156N

## PIEZO SCALER TIPS



88002

Reorder No.155-03  
Endo Chuck



88003

Reorder No.155-01  
Universal A



88004

Reorder No.155-51  
Universal B



88005

Reorder No.155-02  
Chisel



88006

Reorder No.155-48  
Perio

## SPECIFICATIONS

INPUT: AC 100V~240V 50/60Hz

OUTPUT: DC 24V 1.5A

FREQUENCY: 27,000 - 32,000 Hz

DIMENSIONS:  
W221mm x D142mm x H43mm

WEIGHT: 620g

DIMENSIONS(HP): Ø20mmxL 110mm

WEIGHT(HP): 45g

ISO 13485 **GMP**

### DENTAMERICA INC.

18688 E. San Jose Ave. City of Industry, CA 91748 U.S.A.  
TEL: (626) 912-1388 FAX: (626) 913-0510  
[www.dentamerica.com](http://www.dentamerica.com)  
[www.dentamerica.asia](http://www.dentamerica.asia)

DISTRIBUTOR

**DENTAMERICA®**