





EXS[®] Dental Unit

The EXS[®] Dental Unit includes

- Foot Switch for speed control and water activation.
- 3-way air water syringe (air, water and mist).
- High speed handpiece with advanced LED illumination.
- Low speed air motor handpiece system with 1:1 straight attachment for restorative dental attachments and contra angle attachment for prophylaxis.
- Built-in LED Piezo scaler with accessory kit.
- High power suction system.
- Water bottle assembly.
- Quiet and low maintenance oil-free air compressor.



Order No. 102VET(110V) /102DVET(220V)



Easy-to-clean silicone mat and magnetic bur holder keep your workspace tidy and prevent accidental drops.



Intuitive Operation with clearly labeled easy-to-use controls.



Included

Description	Qty
High Speed LED HP	1 pcs
Low Speed Attachment	2 pcs
3-way Air/Water Syringe	1 pcs
Suction System	1 pcs
Disposable Saliva Ejector	1 set (20pcs)
LED Piezo Scaler	1 set
Scaler Tips	3 pcs (Universal*2+Chisel*1)
Prophy Cup/Brush (Latch)	12 pcs each type
Abrasive Polishing Agent	1 btl
High Speed Bur	30 pcs



LED Piezo Scaler Handpiece

Order No. 173H

Straight Attachment

Order No. 166S

Contra Angle Attachment

Order No. 166CA

Low Speed Air Motor Connector

Order No. Z169M

High Speed LED Handpiece

Order No. Z235

High Speed LED Handpiece Connector

Order No. Z236

SPECIFICATIONS of UNIT

Input	AC 100-240V 50/60Hz
Output power consumption	600W
Size	L450 x W600 x H820~920mm
Weight	42KG
Suction Bottle	1500ml

Water Bottle	1000ml
Ultrasonic Scaler Frequency	27,000 to 32,000 Hz
Low speed Attachment speed (Straight, Contra Angle)	20,000 r.p.m.
High speed LED HP speed	380,000 r.p.m.
Compressor Tank Capacity	8L

Apixia Corporation Company Address: 7F-1, 128, Sec.2, Dunhua S. Rd., Da-An Dist., Taipei City Taiwan R.O.C.

Office: +886-2-2702-2608 Fax: +886-2-27554919 Email: vet@apixia.asia www.apixiavet.com

Apixia Inc. Company Address: 18688 San Jose Ave. Industry, CA 91748 USA

Office: +1-626-810-9201 Fax: +1-626-810-9207 Email: info@apixia.com www.apixiavet.com

Apixia Japan Company Address: 1324, Ooida, Daian-cho, Inabe-shi, Mie-ken 511-0274, Japan

Office: +81-0594-41-3080 Fax: +81-0594-41-3081 Email: vet@apixia.asia www.apixiavet.com

Apixia Europe Company Address: 1324, Ooida, Daian-cho, Inabe-shi, Mie-ken 511-0274, Japan

Office: +32-5831-7464 Mobile: +32-4739-37167 Email: pascal@apixia.com www.apixiavet.com