

Dental Hand Instrument 16pcs/set

Order No. 390VET

Winged Elevator

Engineered for precision and ease, our lightweight yet durable winged dental elevators are an essential tool for veterinary dentistry. The innovative winged design enhances grip and control, offering increased leverage for smoother and more efficient tooth extractions, ensuring optimal results for both practitioner and patient.



1mm
Order No.
39001VET



4mm
Order No.
39004VET



2mm
Order No.
39002VET



5mm
Order No.
39005VET



3mm
Order No.
39003VET

Straight Root Tip Pick Stubby Handle

Order No. 39006VET



The straight root tip pick with a stubby handle is designed for optimal grip and control, aiding in the extraction of retained root tips or fractured root fragments during dental surgery.

Tartar Removing Forceps

Order No. 39007VET

These forceps feature a curved beak for easy removal of heavy calculus deposits from the tooth surface prior to scaling, making them an ideal tool for pre-scaling procedures.

Extraction Forceps

Order No. 39008VET

The compact size and spring-loaded handle minimize the risk of applying excessive force, reducing the chance of tooth fractures. These forceps are ideal for veterinary dental procedures, offering both reliability and ease of use.



Minnesota Retractor

Order No. 39009VET

A versatile surgical retractor used to lift and protect soft tissues during high-speed drilling. One end retracts the cheek, while the other holds the tongue away from the surgical field, enhancing safety and visibility during extractions and other oral procedures



Double Ended Periosteal Elevator

Order No. 39010VET

The double-ended periosteal elevator is designed for precise elevation of the mucoperiosteal flap during oral surgery. Its fine, rounded tips are ideal for delicate procedures in cats and dogs and for removing soft tissue or debris from extraction sockets.



Sickle Scaler

Order No. 39011VET

Ideal for removing supragingival plaque and calculus, this scaler features a sharp sickle-shaped tip that is effective in accessing developmental grooves and blood grooves, especially on canine teeth.



Subgingival Curette

Order No. 39012VET

This double-ended curette is used to remove plaque and calculus from the root surface below the gumline. The rounded non-cutting edge minimizes trauma to the gingival tissue while allowing effective subgingival debridement.



Jaquette Scaler

Order No. 39013VET

A double-ended hand scaler designed for removing plaque and calculus above the gumline. Its compact working ends make it particularly effective in tight occlusal spaces and interproximal areas.



Explorer Probe

Order No. 39014VET

Designed for accurate periodontal probing in cats and dogs, featuring 15mm durable etched markings that resist wear from repeated autoclaving.

Luxating Elevator

Specifically crafted for cutting the periodontal ligament and expanding the alveolus, they are perfect for extracting incisors and crowded teeth. The stubby handle ensures safer usage by allowing finger placement near the instrument tip. Two sizes are available for canine and feline teeth.

NOTE: Do not use luxating elevators as winged elevators, and vice versa.



1mm
Order No.
39015VET



2mm
Order No.
39016VET

*Each instrument also available individually.

apixia Corporation Company Address: 7F-1, 128, Sec.2, Dunhua S. Rd., Da-An Dist., Taipei City Taiwan R.O.C.
Office: +886-2-2702-2608 Fax: +886-2-27554919 Email: vet@apixia.asia www.apixiavet.com

apixia Inc. Company Address: 18688 San Jose Ave. Industry, CA 91748 USA
Office: +1-626-810-9201 Fax: +1-626-810-9207 Email: info@apixia.com www.apixiavet.com

apixia Japan Company Address: 1324, Ooida, Daian-cho, Inabe-shi, Mie-ken 511-0274, Japan
Office: +81-0594-41-3080 Fax: +81-0594-41-3081 Email: vet@apixia.asia www.apixiavet.com

BC390VET
Rev.A-Jul/31/2026