

apixia[®] Portable Dental X-ray Unit for VET

Fast, Precise, Low Radiation —
A New Standard in Animal Imaging



Order No. 110VET



Superior Imaging:

0.4mm ultra-fine focal spot for high-contrast details.



Rapid Exposure:

Exposure as fast as 0.04s.



Stable Performance:

High-frequency closed-loop control.



Lightweight & Portable:

Compact and easy to move between rooms.



Intuitive Display:

Clear digital screen shows parameters.






Long Battery Life:

300+ shots per charge.



The apixia® Portable Dental X-ray Unit for VET combines high-frequency technology with an ultra-lightweight design, delivering high-quality images with ease across dental procedures.



ACCESSORIES



Included



Optional

Cone Extension

- Converts the standard cone into a 200mm long-cone configuration for improved positioning consistency.

Backscatter Shield

- 0.53mm lead equivalent
- Helps reduce backscatter radiation for improved operator protection.

SPECIFICATIONS

Model	DIGIREX
Tube Voltage	70kV $\pm 10\%$
Charger Input	100–240V, 50–60Hz
Tube Current	2mA $\pm 20\%$
Focal Spot	0.4mm
Anode Angle	12.5°
Exposure Range	0.04s–2.0s
Relative Humidity	30%–75%RH
Weight	1.9kg

Apixia Corporation Company Address: 7F-1, 128, Sec.2, Dunhua S. Rd., Da-An Dist., Taipei City Taiwan R.O.C.

Office: +886-2-2702-2608 Fax: +886-2-27554919 Email: vet@apixia.asia www.apixiavet.com

Apixia Inc. Company Address: 18688 San Jose Ave. Industry, CA 91748 USA

Office: +1-626-810-9201 Fax: +1-626-810-9207 Email: info@apixia.com www.apixiavet.com

Apixia Japan Company Address: 1324, Ooida, Daian-cho, Inabe-shi, Mie-ken 511-0274, Japan

Office: +81-0594-41-3080 Fax: +81-0594-41-3081 Email: vet@apixia.asia www.apixiavet.com

Apixia Europe

Office: +32-5831-7464 Mobile: +32-4739-37167 Email: pascal@apixia.com www.apixiavet.com